


# Multi-color Panel Design and Experiment Considerations

Ryan Chu 朱伯逢, Ph.D. (c)

Field Applications Scientist

 The world leader in serving science



# Content

1

Basic Principle of Flow Cytometry

2

Attune NxT Flow Cytometer

3

Panel Builder & Experiment Considerations

4

Analyze Extracellular Vesicles by  
Flow Cytometer

5

Benefits and Challenges of Small Particle  
Analysis

6


Tips for Success





# Basic Principle of Flow Cytometry

# Principles of Flow Cytometry




1. Cells in a single profile pass through the narrow tube called flow cell
  - **Cell Focusing**
2. Laser hits individual cell passing through the flow cell
  - **Interrogation Point**
3. Deflected light hits a series of detectors
  - **PMTs**
4. The signals from detectors are interpreted by a computer
  - **Storage and analysis**

What Happens to Light When it Hits a Cell?

## Laser Light Scatter



- When laser light interacts with a cell, light is scattered in all directions
- The magnitude of the light scatter is dependent on refractive index, size and complexity of the cells or particles passing by the laser
- Differences in Forward Light Scatter and Side Light Scatter can be used to distinguish different types of cells or particles



## What Happens to Light When it Hits a Cell?



### Laser Light Scatter

- **Forward Scattered light (FSC)** is impacted by both refractive index and can sometimes be used as a measure of relative cell size.
- **Side-scattered light (SSC)** is a measure of cellular complexity, both surface and internal. SSC is usually collected at 90 degrees to the laser beam.






# Principles of Flow Cytometry




Forward scatter

→ Size



Side scatter


→ Complexity



Fluorescence

→ Phenotype

# What is a Flow Cytometry Antibody?



Conjugated antibodies allow for less experimental steps like washes and centrifugations

Multiple antibodies and reagents are often required for an experiment to detect the cell population of interest



**Fluorophore is attached to the primary antibody**



# Different Fluorophores


- **ORGANIC:**

- Contain several aromatic groups
- FITC, PE, Per-CP, APC, eFluors, Alexa Fluor®, Pacific Blue, Cy5, Cy5.5, Cy7




- **Polymer-dyes**

- Super Bright Dyes, Brilliant Violet Dyes™




- **Non-ORGANIC:**

- Semiconductor particles
- Qdots



# Tandem Dyes (Bi-Molecular Dyes)

## Förster Resonance Energy Transfer (FRET)







Tandem dyes shift emission spectra to higher wavelengths.

Examples:


PE-Cy7, PE-Alexa Fluor® 700, PerCP-Cy5.5, PerCP-eFluor710

APC-Cy7, APC-eFluor 780, Super Bright 600, 645, 702 and 780

# Collect Precise Range of the Emitted Light Wavelengths











## Photomultiplier tube (PMT):

PMT convert photons into electrons and amplify them to create a voltage pulse. Often referred as the “detector”



# Sample Presentation: Voltage Pulse



# Attune NxT Flow Cytometer



# Flow Product line – Analyzer and Cell Sorter Highlight



**Invitrogen™ Attune™  
NxT Flow Cytometer**

2014

Attune NxT Flow Cytometer introduced to the market

2016

Green laser introduced

2017

Invitrogen™ Attune™ violet 6-channel option released

2018

Integration of Thermo Scientific™ Orbitor™ RS2 Microplate Mover

2020

21 CFR Part 11 compliance; Invitrogen™ CytKick™ Autosampler and CytKick™ Max Autosampler with Microsoft™ Windows™ 10 software compatibility released



**Invitrogen™ Attune™ CytPix™  
Flow Cytometer**

2021

Introducing the Attune CytPix Flow Cytometer

2023  
Q2

**Introducing Attune SW 6.0.1 enable AIA (Auto Image analysis)**



**Invitrogen™ Bigfoot Spectral Cell Sorter**

2022

**Introducing Bigfoot Spectral Cell Sorter into APJ**

2023

**New SKU 4L+ SPD activated**

# The Attune NxT Flow Cytometer

## Designed for efficiency, speed and accuracy

With acoustic-assisted hydrodynamic focusing, the Attune NxT Flow Cytometer avoids compromise between data quality and higher sample rates by uncoupling cell alignment from sheath flow. Acoustic-assisted hydrodynamic focusing precisely aligns cells using ultrasonic radiation pressure (>2 MHz) to transport particles into the center of the sample stream. This pre-focused stream is then injected into the sheath stream, resulting in a narrow particle stream and uniform laser illumination, regardless of the sample input rate

### Benefits

- Greater reproducibility and consistency in data
- Maintain consistent concentration results across all flow rates
- Process very dilute or concentrated samples while maintaining low coefficient of variations (CVs)

"The ability to run very dilute samples is quite amazing and might be a life saver on many occasions where you have little-to-no sample left."

– J. P. Robinson, PhD  
Purdue University

acquisition speed

**35,000**  
events/sec

**1000 mL/min**  
max sample input flow rate

Max electronic speed


**65,000**  
events/sec

**10x faster**  
than conventional cytometers



**Speed and accuracy**



# Acoustic Focusing and Hydrodynamic Focusing






# Acoustic Focusing

End-on view of capillary



# Comparable Results at All Flow Rates



		Competitor A		Competitor B		Competitor C		Attune NxT	
	Cells	14 $\mu\text{L}/\text{min}$	66 $\mu\text{L}/\text{min}$	12 $\mu\text{L}/\text{min}$	120 $\mu\text{L}/\text{min}$	12 $\mu\text{L}/\text{min}$	60 $\mu\text{L}/\text{min}$	500 $\mu\text{L}/\text{min}$	1000 $\mu\text{L}/\text{min}$
Seconds	10,000	42.8	9	50	5	50	10	1.2	0.6
Minutes	100,000	7.1	1.5	8.3	0.8	8.3	1.7	0.2	0.1
Minutes	1,000,000	71.3	15.0	83.3	8.3	83.3	16.7	2.0	1.0
Hours	10,000,000	11.9	2.5	13.9	1.4	13.9	2.8	0.3	0.2

# The Attune NxT Flow Cytometer

Transformative, clog-resistant system

## Why Customer Care:

- Clogging number 1 reported issue
- Protect precious samples
- Difficult Sample Types
  - Processed tissue samples, solid tumor samples, adherent cells, HeLa cells, brain tissue, cardiomyocytes, milk, others
- Aggregation tends to occur as cell concentration increases
  - No need to concentrate cells with the Attune NxT, just dilute and run faster!




With the acoustic focusing and sample uptake technology, the Attune rarely clogs or encounters other common fluidic problems as seen with other cytometers that require sheath pressure. During its first year of operation with multiple users we only once required servicing for a problem that was readily solved."

Laurie Kennedy  
Flow Core Manager, University of Calgary

**Protect Precious Samples – Analyze Complex Cells**



# Optical Filters



## Attune NxT accessory filter kits available

- Violet Side Scatter Kit- Small Particles
- Fluorescent Protein Filter Kit- Transfection
- No-wash, No-lyse Filter Kit- Blood

# Choosing Which Channel to Use (BRVY)

Excitation Laser	Emission Filter (nm)	Channel	Common Fluorophores	Fluorescent Proteins
Violet-405 nm	440/50	VL1	SB436, BV421, Pacific Blue, eFluor 450, Alexa Fluor 405	Azurite, Cerulean, eBFP, eCFP, mTurquoise, Sirius
	512/25	VL2	Pacific Green, Qdot 525	T-Sapphire, CFP, vGFP
	603/48	VL3	SB600, BV605, Pacific Orange, Qdot 605	
	710/50	VL4	SB702, Qdot 705	
Blue-488 nm	530/30	BL1	Alexa Fluor 488, FITC, GFP, SYTOX Green	eGFP, Emerald, eYFP
	590/40	BL2	PE, PE-Alexa Fluor 610, PE-Texas Red, Qdot 585	
	695/40	BL3	PE-Alexa Fluor 700, PE-Cy5.5, PerCP, PerCP-Cy5.5, Qdot 705	
Yellow-561 nm	585/16	YL1	Alexa Fluor 568, PE, Qdot 565, PI, SYTOX Orange	mOrange, RFP, dTomato
	620/15	YL2	PE-Alexa Fluor 610, PE-Texas Red	mCherry, DsRed, mKate, mStrawberry
	695/40	YL3	PE-Alexa Fluor 700, PE-Cy5.5, Qdot 705	
	780/60	YL4	PE-Alexa Fluor 750, PE-Cy7, Qdot 800	
Red-637 nm	670/14	RL1	APC, Alexa Fluor 647, SYTOX Red, Qdot 655	
	720/30	RL2	Alexa 700, APC-Alexa Fluor 700, Qdot 705	
	780/60	RL3	APC-Cy7, APC-Alexa Fluor 750, Qdot 800	



## **Panel Builder & Experiment Considerations**

# Multicolor Analysis Workflow



1

Target determination

2

Buffer selection

3

Panel design

4


Biological controls

5

Experimental controls

[thermofisher.com/flowpanel](https://thermofisher.com/flowpanel) for Invitrogen™ Flow Cytometry Panel Builder  
[thermofisher.com/5stepsicfc](https://thermofisher.com/5stepsicfc) for additional resources

# Step 1: Target Determination



- Choose which markers will best answer your research questions
- Which antibodies are needed for markers of interest
- Example T cell panel that includes:
  - Phenotyping markers
  - Cytokines
  - Signaling markers




# Step 1: Target Determination

Phenotyping Markers

T cell identification

T cell subsets



Naive vs. memory T cells



# Step 1: Target Determination


## Signaling Marker

CD3 and CD28 Abs  
for TCR engagement




Staining of phospho-CD247 (Y142) following stimulation with functional anti-CD3 and anti-CD28 antibodies

No staining observed in unstimulated or CD3<sup>-</sup> cells



■ Unstimulated



■ TCR stimulation


# Step 1: Target Determination

## Effector Function

Examples of cytokine staining on stimulated cells

**IL-2:** essential for proliferation and function of activated T cells

**IFN gamma:** anti-viral, anti-tumor cytokine with a wide variety of immune effects



# Step 1: Target Determination

## Optimized Multicolor Immunofluorescence Panels (OMIPs)

- Optimized Multicolor Immunophenotyping Panels (OMIPs) are an additional resource for target selection
- Validated panels for various cell types with gating strategies, clones, and fluorochromes

[Home](#) > [Life Sciences](#) > [Cell Analysis](#) > [Flow Cytometry](#) > [Flow Cytometry Learning Center](#) > [Flow Cytometry Resource Library](#) > [Optimized Flow Cytometry Multiplex Panels](#)

## Optimized Flow Cytometry Multiplex Panels

[Flow Cytometry Resource Library](#)


[Flow Cytometry Application Notes](#)

[Flow Cytometry Research Tools](#)

[Flow Cytometry Educational Videos & Webinars](#)

**Optimized Flow Cytometry Multiplex Panels**

[Flow Cytometry Guided Learning](#)



### Validated antibody panels for flow cytometry

Coined in 2010, the term “optimized multicolor immunofluorescence panel”, or OMIP, refers to a thoroughly tested and validated set of antibodies and reagents that can be used together for the multicolor characterization of a specific cell state or response. Published in the journal *Cytometry Part A* (Wiley Online Library), these OMIPs are designed for flow cytometry, but OMIPs may potentially be defined for image cytometry, fluorescence microscopy, and other polychromatic fluorescence-based methods.


Learn more about Optimized Multicolor Immunofluorescence Panels (OMIPs) in [BioProbes 74](#).

**More targets are not necessarily better!**

# Step 2: Buffer Selection

2

## Buffer selection



Choose the appropriate buffers for sample prep and staining

Intracellular staining always requires fixation and permeabilization of cells



Various options are available, so selecting appropriate option is critical for optimal staining





# Step 2: Buffer Selection

## Optimal Protein Detection

### Intracellular Fixation & Permeabilization Set



### Transcription Factor Staining Buffer Set



### IL-2 staining (cytoplasmic target)

- Optimal with Invitrogen™ eBioscience™ Intracellular (IC) Fixation and Permeabilization Set

### Foxp3 staining (nuclear target)

- Optimal with Invitrogen™ eBioscience™ Foxp3 Transcription Factor Staining Buffer Set

If targeting cytoplasmic IC proteins, first option is most appropriate

Some targets are compatible in both buffers

# Step 2: Buffer Selection


## Why Is Buffer Important?

- The fixation / permeabilization kit allows for access to intracellular targets without damaging epitopes of interest
- General guidelines:
  - Invitrogen™ eBioscience™ Intracellular (IC) Fixation and Permeabilization Set will not allow access to nuclear targets
  - Invitrogen™ eBioscience™ Foxp3 Transcription Factor Staining Buffer Set may not work with some cytoplasmic targets
  - Some phospho-flow targets will require methanol permeabilization (eg. pSTATs)
  - Check product data pages or other online resources
- What if you have targets of interest with different buffer requirements?
  - **Divide experiment into two panels**

# Step 3: Panel Design

3

Panel  
design



Understand differences in fluorochrome performance

Understand the biology of your targets of interest, such as expression levels, co-expression, and gating strategy

Use online resources to help build your panel

# Building Your Panel—Using the Panel Builder




A quick and easy-to-use web tool



Allows incorporation of antibodies using your flow cytometer's configuration



Fluorochrome selection built on spectral visualization of all fluorochromes per laser





[thermofisher.com/flowpanel](https://thermofisher.com/flowpanel)



# Flow Cytometry Panel Builder



Simplified panel design with a 5-step panel design strategy



Find out more at: [thermofisher.com/flowpanel](https://thermofisher.com/flowpanel)

For Research Use Only. Not for use in diagnostic procedures. Not for resale.

# Building Your Panel


## Selecting Your Cytometer

- Extensive list of instruments
- Edit instrument settings or build from scratch
- Remembers settings


Select your cytometer from the drop-down list

Cytometer settings can be customized after selection


### Flow Cytometry Panel Builder<sup>beta</sup>



A simplified immunophenotyping panel design experience including instrument configurations and integrated SpectraViewer plots.



This 5 step process includes a comprehensive summary with easy purchase and export options.





While this beta version is fully functional, additional features and enhancements are coming soon.


STEP 1


Your cytometer  
Attune NxT, Violet(6), Blue(2), Yellow(3), Red(3)

Violet 405nm	Blue 488nm	Yellow 561nm	Red 637nm
450/40	530/30	585/16	670/14
525/50	695/40	620/15	720/30
610/20		780/60	780/60
660/20			
710/50			
780/60			

 Edit cytometer settings

 Load an existing panel


 Clear current panel

 Review SpectraViewer

1: Cytometer   2: Antigens   3: Fluorochromes   4: Products   5: Summary

Save

Next step




This first step lets the panel builder know what fluorochromes your flow cytometer can detect.

# Building Your Panel

## Modify Instrument Settings

Enter cytometer information

LASER		FILTER 	
<div>Laser Name Violet</div>	<div>405</div> <div>nm</div>	<div>V1 450/40</div>	
<div>Remove laser</div>		<div>V2 525/50</div>	<div>×</div>
		<div>V3 610/20</div>	<div>×</div>
		<div>V4 660/20</div>	<div>×</div>
		<div>V5 710/50</div>	<div>×</div>

### Manually specify:

- Lasers
- Fluorochrome/channel names
- Bandpass filters



Editing cytometer from the list and manual editing follow similar designs




Configuration parameters on the cytometers can be saved

Ensure that your filter settings match those of your cytometer so that you can detect your chosen fluorophores

# Building Your Panel

## Selecting Your Antibodies and Reagents



## Slider for Protein Abundance

Keep on medium if you do not know the density


**Antigen density along with fluorochrome brightness are used to recommend fluorochromes**



# Building Your Panel

## Antigen Expression Levels: High


- Markers with **High** expression
  - Well-characterized markers
  - Identify broad subsets of cells
  - Expression is usually positive or negative
  - Often used as parent gate
- Assign these to the “dimmer” fluorochromes, or those most affected by spillover



# Building Your Panel

## Antigen Expression Levels: Medium


- Markers with **Medium** expression
  - Well-characterized markers
  - Defined expression patterns
  - Expression may be on a continuum
- Assign these to the next tier of fluorochromes
  - Intermediate brightness
  - Few spillover problems



# Building Your Panel


## Antigen Expression Levels: Low or Undefined

- Markers with **Low or Undefined** expression
  - There are often fewer color choices for these
  - Often the “critical target” in the panel
- Save your brightest fluorochromes for these!
  - Minimize spillover from other channels
- **Intracellular antigen** expression
  - Staining tends to be dimmer
  - Populations may be less discrete
- Expression levels may vary by donor or activation



# Voltage Optimization

- Also known as voltage walk or voltration
- Why do it?
  - Assures the best sensitivity, resolution, and dynamic range of each PMT detector (and fluorophore)
  - Assures linearity
  - If less than optimal, the resolution of dim populations can be lost
  - More than the optimal voltage gives no advantage to resolution
- When to do it:
  - New machine
  - New filter
  - New alignment
  - New PMT or a Laser
  - New dye (possibly)




**Increasing voltages until there is no increase in resolution**

# Spreading Error

Brighter is Not Always Better

- Spreading Error is spreading of cell populations into other channels
- Spreading Effect is a result of the quantity and the energy of the spillover photons
- Factors that cause Spreading Error
  - Spillover of a fluorophore into adjacent detector (Broad fluorophore emission = Higher Spillover)
  - Balance of Antigen density and Dye Brightness (Bright dye on High expressor = High Spillover)
  - Amplification of photons in adjacent detector (Increasing voltage = Higher Spreading)
  - Energy of Light (Low energy of light = High amplification = High Spreading)


## PE into PerCP-Cyanine5.5



Compensation is not responsible for spread effect,  
compensation only makes spreading visible



# Spreading Error Depends on Several Factors



- Higher compensation values tend to correlate with higher spreading error
- High compensation values are linked to a loss of resolution
- Detector efficiency/amplification of photons = Higher spreading
- Low voltages = Can help with dim population resolution
- Avoid using dim or low stain in channels with high spreading error (Avoid Red and Far Red)


# Antibody Titration – What and Why?

## ➤ What is it?

- Serial dilution of the antibody to determine the concentration for optimal performance
- Can save on reagent usage – compare cost per test
- Recommended working concentrations are just a starting point

## ➤ Why do it?


- Determine target saturation
- Improve resolution
- Eliminate or reduce false positives



Stain Index (SI): 
$$\frac{\text{Median (pos)} - \text{Median (neg)}}{2 \times \text{SD (neg)}}$$

- Calculate staining index to determine best resolution
- Complete for each channel being used

- Understand fluorochrome performance in order to optimize core markers in your panel
- Save the brightest fluorochromes for antigens with low, variable, or unknown expression
- Stain index (SI) defines fluorochrome brightness




$$\text{Stain index} = D/W$$

- D = Difference between positive and negative peak medians
- W = Spread of negative peak, 2 x rSD (robust standard deviation)
- SI is more useful for accessing fluorochrome brightness than the signal-to-background ratio, as it takes W into account

# Brightness of Fluorochromes



Stain Index Ranking



Collected using the Invitrogen™ Attune™ NxT Flow Cytometer, with four laser configuration

PBMCs were labeled with anti-human CD4 antibody conjugated to various fluorochromes, and SI was calculated.

	Fluorochrome	Laser	Filter	Stain Index
Brightest	PE	561	585/16	726
	PE-CY7	561	780/60	565
	PE-CY5.5	488	695/40	542
	BD Horizon PE-CF594	561	620/15	491
	PE-eFluor 610	561	620/15	415
	PE-CY5	488	695/40	393
	PE/DAZZLE 594	561	620/15	358
	BD Horizon BB 515	488	530/30	323
	APC	637	670/14	252
	BD Horizon BV421	405	450/40	233
	Alexa Fluor 647	637	670/14	228
	Super Bright 436	405	440/40	216
	eFluor 660	637	670/14	193
Bright	PERCP-eFluor 710	488	695/40	191
	PE-Alexa Fluor 610	561	620/15	169
	PE TX RED	561	620/15	162
	Alexa Fluor 488	488	530/30	143
	PE-Alexa Fluor 700	488	695/40	138
	BD Horizon BV786	405	780/60	123
	APC-eFluor 780	637	780/60	120
	BD Horizon BV711	405	710/50	111
	Super Bright 702	405	710/50	98
	APC-CY7	637	780/60	95
	Super Bright 780	405	780/60	93
	PERCP-CY5.5	488	695/40	90
	FITC	488	530/30	88
Moderate	Super Bright 600	405	610/20	82
	BD Horizon V450	405	440/40	72
	BD Horizon BV785	405	780/60	72
	APC-Alexa Fluor 750	637	780/60	67
	eFluor 450	405	440/40	67
	Alexa Fluor 405	405	440/40	63
	BD Horizon BV650	405	660/20	63
	BD Horizon BV510	405	512/50	62
	Alexa Fluor 700	637	670/14	54
	Super Bright 645	405	660/20	45
	BD Horizon BV605	405	610/20	42
	APC-H7	637	780/60	35
	Pacific Blue	405	450/40	29
Dim	BD Horizon V500	405	525/50	29
	PERCP	488	695/40	28
	eFluor 506	405	525/50	19
	BD Horizon BV570	405	610/20	15
	Pacific Orange	405	610/20	7

# Building Your Panel




## Selecting Your Fluorochromes

Matrix view of options

Channels: columns

Antigens: rows






# Building Your Panel



## Selecting Your Fluorochromes



Many channels will have more than one fluorochrome option (multiple dots)

### Colored dots

Indicate available fluorochromes in each channel

Black flag: recommended channel


### Recommendations based on:

- Protein abundance from antigen slider
- Fluorochrome brightness

# Building Your Panel



## Filling in Your Panel



### As panel is filled in:

- Colored box indicates chosen channel for the target
- Grayed out column indicates channel is used
- Grayed out row indicates target selection
- Emission curves are filled into SpectraViewer tool as fluorochromes are selected (top)
- Use spillover information from spectraviewer with recommendations based on protein abundance to make the best choices


# Building Your Panel

## Selecting Your Fluorochromes

Tutorial to explain features

Bandpass wavelength

Fluorophore choices



### Black Flag:

Chose “low” protein abundance → bright and medium-bright fluors are flagged

Chose “high” protein abundance → dim and medium-dim fluors are flagged

Chose “medium” or makes no choice → no fluors are flagged


# Building Your Panel

## Selecting Your Fluorochromes

### Black Flags

Flags indicate preferred fluorophores

Brightness indicator for fluorophore



# Building Your Panel




## Visualize Your Fluorochrome Spectrums

### Full-page SpectraViewer on top of page

- Separated by laser
- Spillover values in table under each emission for a given laser line
- Populates as selections are made

Expands to full-screen





# Building Your Panel



## Examine Your Options

Choose the products that work best for you

Select products and packaging size

• CD3, eFluor 450

PRODUCTS (4)	CLONE	TARGET SPECIES
Invitrogen CD3 Monoclonal Antibody (UCHT1), eFluor 450, eBioscience™	UCHT1	Human
Invitrogen CD3 Monoclonal Antibody (OKT3), eFluor 450, eBioscience™	OKT3	Human
Invitrogen CD3 Monoclonal Antibody (SK7), eFluor 450, eBioscience™	SK7	Human Chimpanzee
Invitrogen CD3 Monoclonal Antibody (UCHT1), eFluor 450, eBioscience™	UCHT1	Human

- 1 Antigen/fluorochrome combinations may have multiple product options
- 2 Different clones
- 3 Different sizes
- 4 Clicking on product name takes you to the product data page that includes:

○ Data

○ References

○ Detailed product information

Invitrogen

CD3 Monoclonal Antibody (OKT3), eFluor 450, eBioscience™

✓ Advanced Verification

This Antibody was verified by Relative expression to ensure that the antibody binds to the antigen stated. [View Details](#)

📄 22 Published Figures

📖 24 References

[View \(211\) other CD3 antibodies](#)

# Building Your Panel



## Select Your Clone

Invitrogen

CD3 Monoclonal Antibody (OKT3), Super Bright 600, eBioscience™

28 figures ▼

OKT3

Human

USD 335.00  
Cat # 63-0037-42

Select

100 Tests

16 References

Antibody Testing Data (1)

Advanced Verification (1)

Published Figures (26)

Side Scatter

Forward Scatter

Granulocytes

Monocytes

Lymphocytes

% of Max

CD3 (clone OKT3)

CD3 (clone OKT3)

CD19

80.0

0.48

8.53

11.0

FIGURE 2/28

✓ Relative expression

This Antibody was verified by Relative expression to ensure that the antibody binds to the antigen stated.

CD3 Antibody (63-0037-42)

Staining of human peripheral blood cells. As expected based on known relative expression patterns, CD3 clone OKT3 stains a subset of lymphocytes (T cells) and does not stain monocytes and granulocytes (middle plot). Additional analysis of lymphocytes shows that CD3 clone OKT3 does not stain any CD19+ B cells (right plot). Details: Normal human whole blood was surface stained with CD3 (clone OKT3) and CD19 (... [View more](#))

Internal testing data

Published figures

# Building Your Panel



Export Your List of Antibodies

Order

Export as spreadsheet or PDF

Final panel includes:

- Selected products
- Placeholder reagents
- Other reagents (beads, buffers, and isotype controls)
- SpectraViewer

Can go back to other steps

ATTUNE NXT, VIOLET(6), BLUE(2), YELLOW(3), RED(3)

My Panel

Edit panel

Total: 7 items  
(USD) 2,243.00

Check your price

Add all to cart

.csv

Export as SpreadsheetDownload PDF for printing

Violet Laser

405nm

View SpectraViewer

Channel	Fluorochrome	Product	Price (USD)	Quantity	Select
450/40	eFluor 450	Invitrogen CD4 Monoclonal Antibody (RPA-T4), eFluor 450, eBioscience™	USD 258.00 Cat # 48-0049-42 100 Tests	1	<input checked="" type="checkbox"/>
525/50	LIVE/DEAD Fixable	Invitrogen	USD 281.00	1	<input checked="" type="checkbox"/>

Review SpectraViewer


1: Cytometer2: Antigens3: Fluorochromes4: Products5: Summary

Save

# Step 4: Biological Controls

4

## Biological controls



Biological controls assist with experimental analysis

Help ensure that gates are drawn appropriately


Eliminate false positives and ensure antigen specificity

# Step 4: Biological Controls


## Stimulation Controls

Use unstimulated cells as the control for any cell treatment

Stimulated



Unstimulated




An internal negative control may be more appropriate in some experiments.

**The more controls, the better!**



# Step 4: Biological Controls

## Viability and Singlet Gating




- Dead cells, debris, and coincident events can result in false positives
- Applying a singlet gate removes false positives caused by coincident events
- Using a viability dye excludes dead cells and debris

# Doublet Discrimination

## Examples of Determination


In most cases, data analysis should include gating on single cells.



# Step 4: Biological Controls


## Viability Dye Options

- Exclusion stain, so negative (or dim) population has cells of interest
- Viability dyes detect integrity of the plasma membrane
- Nucleic acid dyes:
  - Can cross compromised membrane to bind nucleic acid
  - Not compatible with cell fixation
- Amine-reactive dyes:
  - Bind surface amines in viable cells (dim staining)
  - Bind IC amines when membrane is compromised (bright staining)
  - Staining is preserved after fixation
- Many options across multiple laser lines
- Can also be part of a dump channel



Live Cells

Dead Cells




LIVE/DEAD™ Fixable Aqua

# Impermeant Nucleic Acid Dyes, Flow Cytometry

Dyes which penetrate cells with a compromised cell membrane to stain nucleic acids, but do not cross the membranes of live cells


- Can be used to **identify dead cells** in a population
- Can be used to **quantitate DNA content** in fixed cells
  - Propidium Iodide (488 nm ex)
  - 7-AAD (488 nm ex)
  - DAPI
  - SYTOX™ AADvanced™ dead cell stain (488 nm ex)
  - SYTOX™ Green dead cell stain (488 nm ex)
  - SYTOX™ Orange dead cell stain (488/532/561 ex)
  - SYTOX™ Blue dead cell stain (405 nm ex)
  - SYTOX™ Red dead cell stain (633 nm ex)



# SYTOX™ Dead Cell Stains



Five different colors for flexibility in multicolor panels



# Propidium Iodide ReadyProbes™ Reagent

Power and Simplicity




- Ready-to-use liquid propidium iodide formulation
- Rapid staining of dead cells without wash steps
- Ready-to-use liquid formulation in convenient dropper bottle—no need to dilute, weigh, or pipette
- **Stable at room temperature**—keep handy at your work station or cell culture area

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.  
For Research Use Only. Not for use in diagnostic procedures. © 2013 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. MAN0077238, Rev. 1.0




# LIVE/DEAD™ Fixable Dead Cell Stains



## Live cells:

React with the fluorescent reactive dye only on their surface to yield weakly fluorescent cells.




## Cells with compromised membranes:


The dye reacts with free amines both in the cell interior and on the cell surface, yielding intense fluorescent staining.

**Viability** = membrane integrity

## Before fixation



## After fixation



Fixable violet dead cell stain  
405 nm Violet Excitation (440/40 BP)

# LIVE/DEAD™ Fixable Dead Cell Stains

- **Compatible with Fix and Perm Procedures**

Useful for dead cell identification with intracellular targets

- **Specificity and reliability**

Staining pattern is similar before and after fixation


- **Reduced handling risks**

Formaldehydes are known to reduce risks of viruses


- **Convenience**

Cells can be stained and fixed at various times during the experiment, and the results can be analyzed several hours later, without loss of the discrimination pattern

## Propidium Iodide





## LIVE/DEAD™ Fixable Red Stain



# Amine Reactive Dyes for Dead Cell Identification



8 Color Options, Packaged for Stability






## LIVE/DEAD™ Fixable Dead Cell Stain Kits

Aqua-fluorescent reactive dye

Reactive dye	Excitation source	Ex*	Em*
Blue fluorescent reactive dye (L23105)	UV	350	450
Violet fluorescent reactive dye (L34955)	405 nm	416	451
Aqua fluorescent reactive dye (L34957)	405 nm	367	526
Yellow fluorescent reactive dye (L34959)	405 nm	400	575
Green fluorescent reactive dye (L23101)	488 nm	495	520
Red fluorescent reactive dye (L23102)	488 nm	595	615
Far red fluorescent reactive dye (L10120)	633/635 nm	650	665
Near-IR fluorescent reactive dye (L10119)	633/635 nm	750	775

\* Approximate fluorescence excitation (Ex) and emission (Em) maxima, in nm.



## Impermeant DNA dyes

- Add at final step, do not wash out
- Emission is broad, consider for multicolor applications
- Dead/Fixed cells can be used for compensation control


## Amine-reactive dyes

- Do not use protein in buffers
- Live cells have dim fluorescence
- Use with -aldehyde fixatives
- Can be used without fixing cells too
- ArC™ compensation beads useful

# Step 5: Experimental Controls


- **Unstained cells:** to assess AF of cells
- **Single-color controls:** to set compensation and optimize voltages
- **FMO controls:** to assess the effect of other fluorochromes and compensation on background
- **Isotype controls:** to assess the stickiness of antibodies
- **Biological controls:** viability for removing dead cells and debris; unstimulated for correct gating

The more controls, the better!




— Unstained — FMO — Isotype control — Experimental Ab


# Compensation: Eliminates Spectral Overlap



Un-Compensated  
Single Color Control



Compensated  
Single Color Control



**R-PE channel requires more compensation than Alexa Fluor 488 channel**



# Basic Rules of Compensation

- Unstained cells
- Single color controls are required
- Controls need to be at least as bright as the brightest positive sample
- Background fluorescence should be the same for the positive and negative control populations
- Compensation color must be matched to your experimental color (FITC cannot substitute for GFP)
- The actual tandem dye being used in the sample staining must be used in the single-color control
- Collect enough events to be statistically relevant


# Compensation Beads

Do I need compensation beads?

- Cells may be in limited supply – reserve for experimental use
- Antigen may be dimly expressed by the cells
- Compensate with the antibody used in experiment
- Few cells in sample express the antigen of interest
- Correct compensation guaranteed with positive & negative populations available

# Are You Simply Guessing?

Good compensation controls contain distinct positive and negative populations



IL-13 staining on cells

- Broad distribution
- Few positive events
- Difficult to identify positive population

IL-13 staining on OneComp eBeads

- Discrete distribution
- 50% positive
- Easy to identify positive population

# Compensation Beads



	UltraComp eBeads Plus compensation beads	UltraComp eBeads compensation beads	OneComp eBeads compensation beads	AbC™ Total Antibody Compensation Bead Kit
Reactivity	Human, rabbit, hamster, mouse, and rat antibodies with recognition of the kappa and lambda chains.	Hamster, mouse, and rat antibodies with recognition of the kappa and lambda chains.	Hamster, mouse, and rat antibodies with recognition of the kappa and lambda chains.	Mouse, rat, hamster, and rabbit antibodies with recognition of the kappa and lambda chains.
Laser compatibility	UV to 633 nm Improved for polymer dye use from violet laser	UV to 633 nm	Not with UV or violet lasers.	UV to 633 nm

OneComp eBeads Compensation Simplified

One drop does it all


One Drop Does it All

UltraComp eBeads™


The Ultimate Compensation Bead




# Fluorescence Minus One (FMO)



- Setting gates is not arbitrary!
- Don't rely only on unstained
- Use FMO to take into account spreading from compensation
- Set up FMO for every color for pilot study
- Especially important when +/- populations are not easily resolved



# Un-Stimulated Controls Inform on Gating Strategy



Unstimulated cells is an important biological control for gating strategy



## Steps:

1

### Target determination

Choose the phenotyping and functional markers to analyze your cells of interest

2

### Buffer selection

Determine which buffers will provide optimal staining

3

### Panel design

Select fluorochromes to optimize core markers and minimize spillover

4

### Biological controls

Include viability and stimulation controls to remove dead cells and analyze accurately

5


### Experimental controls

Set up experimental controls to compensate for spillover

A woman with dark hair, wearing a light blue surgical face mask and a red and white patterned scarf, is looking down at a smartphone in her right hand. She is wearing a grey jacket. The background is a blurred city street at night with warm bokeh lights from street lamps and other pedestrians. The overall tone is cool and blue.

## Analyze Extracellular Vesicles by Flow Cytometer

# Flow cytometry: what does it mean?



# CYTOMETRY



Cell



Measurement

measurement of a ~~cell~~ particle

# EV formation and release

## Extracellular Vesicles

- Family of small membrane vesicles
- Found in the extracellular environment
- Carry information about their cell of origin
- Phospholipid bilayer
- Found in all body fluids
- Found in cell cultures

## Exosomes

- Multi-Vesicular Body (MVB) formation
- MVB outer membrane fuse with the cell membrane
- Vesicles inside the MVB exit the cell
- Outside of the cell membrane is on the outside of the exosome

## Extracellular Vesicles

```
graph TD; EV[Extracellular Vesicles] --> EX[Exosomes]; EV --> MV[Microvesicles];
```

## Microvesicles

- Outward budding of the cell membrane
- Outside of the cell membrane is on the outside of the microvesicle

## Apoptotic Bodies

- Cell fragmentation, often called blebbing

Journal of Extracellular Vesicles

*Minimal information for studies of extracellular vesicles 2018 (MISEV2018): a position statement of the International Society of Extracellular Vesicles*

<https://doi.org/10.1080/20013078.2018.1535750> Open access

Guidelines for Flow Cytometry

*MIFlowCyt-EV: a framework for standardized reporting of extracellular vesicle flow cytometry experiments*

➤ ISEV, ISAC and ISTH joint working group EV guidelines


<https://www.tandfonline.com/doi/full/10.1080/20013078.2020.1713526> Open access

**“ISEV endorses “extracellular vesicle” (EV) as the generic term for particles naturally released from the cell that are delimited by a lipid bilayer and cannot replicate.**

Since consensus has not yet emerged on specific markers of EV subtypes, such as endosome-origin “exosomes” and plasma membrane-derived “ectosomes” (microparticles/microvesicles) assigning an EV to a particular biogenesis pathway remains extraordinarily difficult unless, e.g. the EV is caught in the act of release by live imaging techniques.”

# Extracellular Vesicles increasing interest


- Interest in studying EVs has grown exponentially in the last decade
- EVs come from many different cell types and fluid sources
  - Cultured Cells, Tumor and Cancer Cells, Normal Cells
  - BioFluids (plasma, urine, CSF)
- Evidence that EVs play an important role in health and disease
  - Involved in intercellular communication
  - May transfer cargo
- Measurement of EVs and their cargo help in understanding their origins, potential targets, and possible functions







# Extracellular Vesicles increasing interest

- EVs can be loaded with multiple proteins, nucleic acids, DNA, RNAs (mRNA, miRNA, ncRNA), lipids, and other molecules (i.e., biomolecules of viruses)
- Represent potential diagnostic and prognostic targets
  - Biomarker for cancer and immunosuppression
  - Cardiovascular disease
  - Infectious diseases
  - Inflammation




# Relative sizes



György, B., Szabó, T.G., Pásztói, M. et al. Cell. Mol. Life Sci. (2011) 68: 2667.  
<https://doi.org/10.1007/s00018-011-0689-3>

# Extracellular vesicles



	Exosomes	Microvesicles	Apoptotic Bodies
Size	30-100 nm	100-1000 nm	1000-5000 nm
Shape	regular	irregular	irregular
Origin	Endosomal compartments	Cell Surface Plasma Membrane	Cells undergoing Apoptosis


Adapted from: Oncology Reports 35(35):1237-1224 Dec 2015

A woman with dark hair, wearing a light blue surgical face mask and a red and white patterned scarf, is looking down at a smartphone in her right hand. She is wearing a grey jacket. The background is a blurred city street at night with warm bokeh lights from street lamps and buildings. Other pedestrians are visible in the background, also blurred.

## Benefits and Challenges of Small Particle Analysis



# Extracellular Vesicle characterization methods



Jean-Daniel Tisson, et al, Translational Proteomics 2013, <https://doi.org/10.1016/j.trprot.2013.04.004>

# Extracellular Vesicle characterization methods



Two Categories:

- 1. Total or aggregate EV measurements
- 2. Individual EV measurements: important in understanding heterogeneity

Characterization Method	Assay Type	Measurements
Flow Cytometry	Individual and Total Vesicles	Concentration, Cargo, Size*
Nanoparticle Tracking Analysis (NTA)	Individual Vesicles	Concentration**, Size
Resistive Pulse Sensing (RPS)	Individual Vesicles	Concentration**, Size
Dynamic Light Scattering (DLS)	Individual Vesicles	Size
Atomic Force Microscopy	Individual Vesicles	Visual Confirmation, Size
Electron Microscopy (EM)	Individual Vesicles	Visual Confirmation, Size
Western Blot	Total Vesicles	Cargo
ELISA	Total Vesicles	Cargo
* with appropriate size standards		
**Particles, not vesicle specific		



# Benefits of flow cytometry

- Records data from single particles
- Makes measurements on large numbers of particles
- Thousands of particles can be analyzed rapidly
- Rich statistical analysis on populations
- Because single particles are measured, heterogeneity can be revealed
- Ability to archive standardized files
- Flexibility of data analysis with the ability to re-analyze data
- Ability for automated high throughput acquisition
- Ability to multiplex
- Counting ability

# Benefits of flow cytometry in small particle testing

## Single particle analysis

- Ability to analyze small particle rapidly and individually
  - Traditional methods of EV analysis are limited in the capacity to analyze individual particles
  - Time-consuming and expensive options limit EV studies
- Capable of analyzing thousands of particles per second, improving statistics, and allowing the quantification of unique and rare particle types.

## Particle counting and concentration

- A significant advantage in both speed and accuracy for counting methods
- Obtain counts by direct volumetric counting with the Invitrogen™ Attune™ NxT Flow Cytometer

## Ability to Multiplex

- More information can be obtained by simultaneous analysis with fluorescent labels such as surface proteins, membrane lipids, cellular esterase


# Challenges with EV detection using flow cytometry

## Identity



Small size of EVs, their heterogeneity, and low refractive index make them difficult to identify<sup>1,3</sup>

## High Coincidence



Multiple small particles can be erroneously detected as single microvesicles<sup>1,2</sup>

## Detection



Fluorescent labeling of EVs required for optimal detection<sup>1</sup>


1. Nolan, J.P. 2015. Flow cytometry of extracellular vesicles: potential, pitfalls, and prospects. *Curr. Protoc. Cytom.* 73:13.14.1-13.14.16. doi: 10.1002/0471142956.cy1314s73
2. Rousseau M, Belleannee C, Duchez A-C, et al. Detection and Quantification of Microparticles from Different Cellular Lineages Using Flow Cytometry. Evaluation of the Impact of Secreted Phospholipase A<sub>2</sub> on Microparticle Assessment. Combes V, ed. *PLoS ONE*. 2015;10(1):e0116812.
3. Pospichalova V, Svoboda J, Dave Z, et al. Simplified protocol for flow cytometry analysis of fluorescently labeled exosomes and microvesicles using dedicated flow cytometer. *Journal of Extracellular Vesicles*. 2015;4:10.3402/jev.v4.25530.

## Tips for Success

# Anatomy of a light pulse

Particles flow through a focused spot of light

Sample Stream




The diagram illustrates the interaction between a sample stream and a light pulse. On the left, a vertical black arrow points downwards, labeled 'Sample Stream'. A small red dot is positioned at the top of this arrow. In the center of the diagram, there is a larger red dot. Directly beneath this red dot is a horizontal, glowing green oval, representing a focused spot of light. To the right of the green oval is a large, solid white rectangle, which likely represents a detector or a collection area. The entire scene is set against a light gray background.

This generates a pulse of light

# Anatomy of a light pulse

The Gaussian-shaped light pulses are measured



Data is derived from these light pulses


The Area, Height, and Width refer to the optical pulses




# Threshold in flow cytometry

Threshold is an electronic hurdle that establishes a criteria whether events are recorded or not

**Above threshold:  
A pulse is recorded**



**Below threshold:  
pulse Not recorded**



# Threshold used in small particle analysis

If every single particle passing through the laser caused the instrument to collect data, the data would be dominated by information coming from a large number of minute particles.

To prevent this, a threshold is set such that a certain size or fluorescent intensity must be exceeded for the instrument to collect data.

Thresholds may be set on a scatter parameter, or a fluorescence parameter, or a combination of parameters.

Threshold

OR FSC 25 x1000


OFF x1000

OFF x1000


OFF x1000

☐ Exclude coincident events


# Forward scatter threshold




# Forward scatter threshold



# Forward scatter threshold




# Forward scatter threshold





# Fluorescence threshold


Increasing Fluorescence Threshold  
Invitrogen™ SYTO™ 9 Green Fluorescent Nucleic Acid stain (BL1-H)




Adjustment of fluorescence threshold is done empirically

- While acquiring sample, and
- Before Recording sample

# Linear vs Log Display for FSC and SSC



Linear scale



Use Log Display  
for small particle  
detection

# Pre-analytic considerations

Pre-analytic steps to consider:

- Source material
  - Collection
  - Transport
  - Storage
- Isolation/enrichment technique
- Exosome-depleted products
- Filter fluids
  - Sheath Fluid
  - Buffers
- Choice of Ab/fluorophore/dye
  - Size
  - Brightness
  - Specificity
- Controls

## Sample Preparation



Journal of Immunological Methods

Volume 438, November 2016, Pages 11–20



Research paper

The impact of various preanalytical treatments on the phenotype of small extracellular vesicles in blood analyzed by protein microarray

Rikke Bæk  , Evo K.L. Søndergaard, Kim Varming, Malene M. Jørgensen

Department of Clinical Immunology, Part of Extracellular Vesicle Research Center Denmark (EVsearch.dk), Aalborg University Hospital, Aalborg, Denmark

Received 5 April 2016, Revised 23 August 2016, Accepted 23 August 2016, Available online 24 August 2016



Gibco™ Fetal Bovine Serum, exosome-depleted

# Sample Isolation methods

## EV Isolation, aka Separation, Purification, Enrichment, Concentration

- Usually required for subsequent analysis
- Each method has its strengths and limitations

## Primary Isolation Methods

- Ultra-centrifugation
- Combinations of Methods
- Density Gradient
- Filtration
- Size Exclusion Chromatography
- Precipitation
- Bead Capture



Total Exosome Isolation Reagent



Journal of  
Extracellular Vesicles

ORIGINAL RESEARCH ARTICLE

## Techniques used for the isolation and characterization of extracellular vesicles: results of a worldwide survey



Chris Gardiner<sup>1\*</sup>, Dolores Di Vizio<sup>2</sup>, Susmita Sahoo<sup>3</sup>, Clotilde Théry<sup>4</sup>, Kenneth W. Witwer<sup>5</sup>, Marca Wauben<sup>6</sup> and Andrew F. Hill<sup>7</sup>

<sup>1</sup>Haemostasis Research Unit, Research Department of Haematology, University College London, London, UK; <sup>2</sup>Icahn School of Medicine, Cardiovascular Research Center, Cedars-Sinai, Los Angeles, CA, USA; <sup>3</sup>Cardiovascular Research Center, Icahn School of Medicine at Mount Sinai, New York, NY, USA; <sup>4</sup>Institut Curie, PSL Research University, INSERM U932, Paris, France; <sup>5</sup>Johns Hopkins University School of Medicine, Baltimore, USA; <sup>6</sup>Department of Biochemistry and Cell Biology, Faculty of Veterinary Medicine, Utrecht University, Utrecht, The Netherlands; <sup>7</sup>Department of Biochemistry and Genetics, La Trobe Institute for Molecular Science, La Trobe University, Bundoora, Australia

Chris Gardiner, et. al., Techniques used for the isolation and characterization of extracellular vesicles: results of a worldwide survey, Journal of Extracellular Vesicles (2016) 5:1, DOI: [10.3402/jev.v5.32945](https://doi.org/10.3402/jev.v5.32945)  
[OPEN ACCESS](#)


# Differential Ultra-Centrifugation

From cell culture supernatant or body fluids



# Dilution for single cell detection

- Higher particle concentrations result in coincidence effect (swarm)
- Significant sample dilution is required to ensure single particle detection






# Flow cytometry challenges

EVs have much smaller surface area and volume compared to cells

Several approaches are needed

- Signal Detection
  - Bright Fluorophores
  - More sensitive SSC
- Background Reduction
  - Attention to system cleanliness
  - UltraFilter Instrument Fluids
  - UltraFilter Buffers
- Acquisition/Analysis
  - Coincidence
  - Thresholds



Extracellular Vesicle  
100 nm diameter

Surface marker density  
~100 per cell

One-millionth volume of  
a 10 µm diameter cell

- Labeling Reagent Selection
  - Bright fluorophore
- Antibodies
  - Abundant targets
  - CD9, CD63, CD81 standard



## Label all EVs (Fluorescence Threshold)

- Lipid specific Di-8-ANNEPS dye
- Lipophilic styryl compounds
  - FM™ 1-43 dye or FM™ 4-64 dye
- Carbocyanine lipophilic dyes
  - Vybrant™ DiO, DiI, DiD
- Invitrogen™ Bodipy™ dyes
  - Bodipy 493/503
  - Bodipy 581/591
- Invitrogen™ CellTrace™ dyes
  - Numerous ex/em

# Additional tips for small particle success

- Background Reduction
  - Use a Fluorescence Threshold
  - Advanced Settings
    - Decrease Windows Extensions
    - Decrease Width Threshold
  - Focusing Fluid with 0.1  $\mu\text{m}$  or smaller Filter
    - ideally filter to 0.03- 0.05  $\mu\text{m}$
  - Filter buffers with 0.1  $\mu\text{m}$  or smaller Filter
  - Change Focusing Fluid Filters Monthly
- Signal Increase



Choice™ PES Syringe Filters



Attune™ NxT Focusing Fluid Filters



Nalgene™ Rapid-Flow™ Sterile Disposable Filter Units with PES Membrane

# Guidelines for acquisition & Concentration/Dilution testing



Example sample type	Sample type	
	0.2–3 µm	>3 µm
	Bacteria Microspheres  Extracellular Vesicles may require additional instrument prep	Jurkat cells Ramos cells Leukocytes Microspheres Cardiomyocytes
Flow rate	12.5–1,000 µL/min	100–1,000 µL/min
Sample concentration	500–10 <sup>6</sup> particles/mL	
Event rate	<8,000 events/sec	
Sample volume	50–4,000 µL	

## Recommendation for EVs:

- 12.5 or 25 µL/min flow rate
- Measure ≥ 100 µL Sample Volume
- Record consistent volume every sample
- Keep Event Rate < 8,000 events/sec
- Attention to cleaning: Instrument & between samples
- Ultra-Filter fluids used, rinse containers

## Ensure Single Particle Analysis-Eliminate Coincidence


Dilution	Theoretical Events/sec	Observed Events/Sec
1:100		5711
1:200	2816	2881
1:400	1408	1568
1:800	704	750

Event Rate decrease in proportion to dilution

- Scatter standardization for sizing of particles relies on theoretical modeling combined with analysis of beads of known diameter and refractive index.
- Availability of reference particles with refractive index close to EVs for standardization would yield accurate EV diameter data directly, without the need for complex models.
  - Silica beads, hollow silica beads, fluorescent silica beads
  - FP labeled EVs
  - Liposomes of known diameter
  - Combination of fluorescent polystyrene beads with silica beads
- Limited commercial availability of standards, but growing

# Refractive Index and scatter

Scatter Intensity




	Liposome	Extracellular Vesicle	Silica Bead	Polystyrene Bead
Refractive Index	1.39	1.37-1.42	1.4	1.59-1.63

NOTE: Use of Polystyrene beads for Scatter ‘size calibration’ is not recommended due to differences in RI. Instead use a particle closer in RI to the EV; helpful to characterize your instrument too.



## Suggested Controls

- Buffer only
- Buffer + all reagents (no EVs)
- Buffer + EVs only (no stain)
- Sample dilutions of EVs to ensure single particle detection
- Particle for instrument characterization
- Optional: isotype controls
- FMO, single color controls
- Detergent treated for negative control ( such as Triton™ X-100 or NP-40)
- Particles for size and fluorescence standardization



# Attune NxT Small Particle Side-Scatter Filter




- Designed to increase the dynamic range of the SSC without sacrificing resolution
- The use of this 488/10 side scatter filter combined with ultra-filtering of the focusing fluid has been found to dramatically decrease background noise and enable discrimination of 100-nm particles
- Do not use the sp-BSSC filter for PT, or for cells
- ID the sp-BSSC in the SSC label
- Signal increase approx. 100-fold



# Attune NxT Small Particle Side-Scatter Filter testing


Ultra-Filtration of Focusing Fluid, running ApogeeMix particles:

- Silica beads sizes: 180, 240, 300, 590, 880, 1300 nm (non-fluorescent)
- ApogeeMix centrifuged and suspended in Ultra-Filtered PBS



Data from Steve McClellan, Mitchell Cancer Institute, University of South Alabama

# ApogeeMix particles with SP Side-Scatter Filter





Silica Particle Sizes (nm)
180
240
300
590
880
1300

RI Silica: 1.43  
RI polystyrene: 1.59

Data from Steve McClellan, Mitchell Cancer Institute, University of South Alabama

# EV testing FM 1-43 dye




Data from Steve McClellan, Mitchell Cancer Institute, University of South Alabama


# EV multi-color testing



Gated on FM 1-43 Positive



Sequential Gate from FM 1-43 Positive,  
then CD63 PE-Cyanine7 Positive



Data from Steve McClellan, Mitchell Cancer Institute, University of South Alabama



# Best practices to detect EVs using flow cytometry

- Attention to pre-analytical steps
- Maintain clean instrument
- Attention to cleaning between samples
- Particles for instrument characterization that have a Refractive Index close to EVs
- Sample flow rate set at low speed
- Fluorescence threshold for EV identification, OR
- Filter all fluids including buffer & sheath fluids
- Sample dilution to minimize coincidence
- Treat samples with detergent to confirm EVs
- Bright fluorophore selection
- Use Attune NxT Small Particle Side-Scatter Filter



# Thank you

[support.tw@thermofisher.com](mailto:support.tw@thermofisher.com)

Rescue Lens App for remote support

